

Product Datasheet

Human LRRK2 ELISA Kit (orb1198810)

Description	Human Leucine Rich Repeat Kinase 2 (LRRK2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q5S007
Sensitivity	0.122ng/mL
Entrez	120892
Expiration Date	Please enquire.