

Product Datasheet

Human LRP4 ELISA Kit (orb1198814)

Description	Human Low Density Lipoprotein Receptor Related Protein 4 (LRP4) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	075096
Sensitivity	0.103ng/mL
Entrez	4038
Expiration Date	Please enquire.