

Product Datasheet

Human LRCH4 ELISA Kit (orb1198817)

Description	Human Leucine-Rich Repeat And Calponin Homology Domain-Containing Protein 4 (LRCH4) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	075427
Sensitivity	27pg/mL
Entrez	231798
Expiration Date	Please enquire.