

Product Datasheet

Human LPCAT3 ELISA Kit (orb1198818)

Description	Human Lysophosphatidylcholine Acyltransferase 3 (LPCAT3) ELISA Kit
Reactivity	Human
Range	3.12-200ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q6P1A2
Sensitivity	1.29ng/mL
Entrez	10162
Expiration Date	Please enquire.