

Product Datasheet

Human LPCAT2 ELISA Kit (orb1198819)

Description	Human Lysophosphatidylcholine Acyltransferase 2 (LPCAT2) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q7L5N7
Sensitivity	0.58ng/mL
Entrez	54947
Expiration Date	Please enquire.