

Product Datasheet

Human LPCAT1 ELISA Kit (orb1198820)

Description	Human Lysophosphatidylcholine Acyltransferase 1 (LPCAT1) ELISA Kit
Reactivity	Human
Range	3.12-200ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q8NF37
Sensitivity	1.15ng/mL
Entrez	79888
Expiration Date	Please enquire.