

Product Datasheet

Human LOXL3 ELISA Kit (orb1198824)

Description	Human Lysyl Oxidase Like Protein 3 (LOXL3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P58215
Sensitivity	0.072ng/mL
Entrez	84695
Expiration Date	Please enquire.