

Product Datasheet

Human LEFTY1 ELISA Kit (orb1198842)

Description	Human Left/Right Determination Factor 1 (LEFTY1) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	075610
Sensitivity	0.24ng/mL
Entrez	10637
Expiration Date	Please enquire.