

Product Datasheet

Human LECT1 ELISA Kit (orb1198844)

Description	Human Leukocyte Cell Derived Chemotaxin 1 (LECT1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	075829
Sensitivity	0.051ng/mL
Entrez	11061
Expiration Date	Please enquire.