

Product Datasheet

Human LDLRAP1 ELISA Kit (orb1198845)

Description	Human Low Density Lipoprotein Receptor Adaptor Protein 1 (LDLRAP1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q5SW96
Sensitivity	0.056ng/mL
Entrez	26119
Expiration Date	Please enquire.