

Product Datasheet

Human ITIH2 ELISA Kit (orb1198899)

Description	Human Inter Alpha-Globulin Inhibitor H2 (ITIH2) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P19823
Sensitivity	35pg/mL
Entrez	3698
Expiration Date	Please enquire.