

Product Datasheet

Human ITIH1 ELISA Kit (orb1198900)

Description	Human Inter Alpha-Globulin Inhibitor H1 (ITIH1) ELISA Kit
Reactivity	Human
Range	31.2-2000ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P19827
Sensitivity	15.7ng/mL
Entrez	3697
Expiration Date	Please enquire.