

Product Datasheet

Human HTR4 ELISA Kit (orb1198924)

Description	Human 5-Hydroxytryptamine Receptor 4 (HTR4) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q13639
Sensitivity	0.063ng/mL
Entrez	3360
Expiration Date	Please enquire.