

Product Datasheet

Human HSPA5 ELISA Kit (orb1198936)

Description	Human Heat Shock 70kDa Protein 5 (HSPA5) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P11021
Sensitivity	0.117ng/mL
Entrez	3309
Expiration Date	Please enquire.