

Product Datasheet

Human HS6ST2 ELISA Kit (orb1198945)

Description	Human Heparan Sulfate-6-O-Sulfotransferase 2 (HS6ST2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q96MM7
Sensitivity	0.066ng/mL
Entrez	90161
Expiration Date	Please enquire.