

Product Datasheet

Human HNRPK ELISA Kit (orb1198956)

Description	Human Heterogeneous Nuclear Ribonucleoprotein K (HNRPK) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P61978
Sensitivity	0.107ng/mL
Entrez	3190
Expiration Date	Please enquire.