

Product Datasheet

Human GUCA2B ELISA Kit (orb1198988)

Description	Human Guanylate Cyclase Activator 2B (GUCA2B) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q16661
Sensitivity	6.2pg/mL
Entrez	2981
Expiration Date	Please enquire.