

Product Datasheet

Human GPAM ELISA Kit (orb1199015)

Description	Human Glycerol-3-Phosphate Acyltransferase, Mitochondrial (GPAM) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9HCL2
Sensitivity	0.059ng/mL
Entrez	57678
Expiration Date	Please enquire.