

Product Datasheet

Human GDF3 ELISA Kit (orb1199041)

Description	Human Growth Differentiation Factor 3 (GDF3) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9NR23
Sensitivity	34pg/mL
Entrez	9573
Expiration Date	Please enquire.