

Product Datasheet

Human FPN ELISA Kit (orb1199069)

Description	Human Ferroportin (FPN) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9NP59
Sensitivity	0.126ng/mL
Entrez	30061
Expiration Date	Please enquire.