

Product Datasheet

Human FLNC ELISA Kit (orb1199077)

Description	Human Filamin C Gamma (FLNC) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q14315
Sensitivity	0.126ng/mL
Entrez	2318
Expiration Date	Please enquire.