

Product Datasheet

Rat ERK1 ELISA Kit (orb1199119)

Description	Rat Extracellular Signal Regulated Kinase 1 (ERK1) ELISA Kit
Reactivity	Rat
Range	0.625-40ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P21708
Sensitivity	0.249ng/mL
Entrez	50689
Expiration Date	Please enquire.