

Product Datasheet

Human ErbB4 ELISA Kit (orb1199121)

Description	Human V-Erb A Erythroblastic Leukemia Viral Oncogene Homolog 4 (ErbB4) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q15303
Sensitivity	0.062ng/mL
Entrez	2066
Expiration Date	Please enquire.