

Product Datasheet

Human ENPP2 ELISA Kit (orb1199131)

Description	Human Ectonucleotide Pyrophosphatase/Phosphodiesterase 2 (ENPP2) ELISA Kit
Reactivity	Human
Range	93.7-6000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q13822
Sensitivity	32.28pg/mL
Entrez	5168
Expiration Date	Please enquire.