

Product Datasheet

Human ECM1 ELISA Kit (orb1199151)

Description	Human Extracellular Matrix Protein 1 (ECM1) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q16610
Sensitivity	13.2pg/mL
Entrez	1893
Expiration Date	Please enquire.