

## Product Datasheet

### Human EAAT2 ELISA Kit (orb1199155)

|                                |                                                                 |
|--------------------------------|-----------------------------------------------------------------|
| <b>Description</b>             | Human Excitatory Amino Acid<br>Transporter 2 (EAAT2) ELISA Kit  |
| <b>Reactivity</b>              | Human                                                           |
| <b>Range</b>                   | 0.312-20ng/mL                                                   |
| <b>Tested<br/>Applications</b> | ELISA                                                           |
| <b>Note</b>                    | For research use only                                           |
| <b>Application notes</b>       | sandwich                                                        |
| <b>Sample Types</b>            | tissue homogenates, cell lysates or other<br>biological fluids. |
| <b>Uniprot ID</b>              | <b>P43004</b>                                                   |
| <b>Sensitivity</b>             | 0.113ng/mL                                                      |
| <b>Entrez</b>                  | <b>6506</b>                                                     |
| <b>Expiration Date</b>         | Please enquire.                                                 |