

Product Datasheet

Human DPAGT1 ELISA Kit (orb1199186)

Description Human Dolichyl Phosphate-N-Acetylglucosaminophosphotransferase 1 (DPAGT1) ELISA Kit

Reactivity Human

Range 0.156-10ng/mL

Tested Applications ELISA

Note For research use only

Application notes sandwich

Sample Types tissue homogenates, cell lysates or other biological fluids.

Uniprot ID [Q9H3H5](#)

Sensitivity 0.059ng/mL

Entrez [1798](#)

Expiration Date Please enquire.