

Product Datasheet

Human DNAH11 ELISA Kit (orb1199194)

Description	Human Dynein, Axonemal, Heavy Chain 11 (DNAH11) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q96DT5
Sensitivity	0.127ng/mL
Entrez	8701
Expiration Date	Please enquire.