

## Product Datasheet

### Human DIAPH1 ELISA Kit (orb1199205)

|                            |                                                |
|----------------------------|------------------------------------------------|
| <b>Description</b>         | Human Diaphanous Homolog 1 (DIAPH1) ELISA Kit  |
| <b>Reactivity</b>          | Human                                          |
| <b>Range</b>               | 0.156-10ng/mL                                  |
| <b>Tested Applications</b> | ELISA                                          |
| <b>Note</b>                | For research use only                          |
| <b>Application notes</b>   | sandwich                                       |
| <b>Sample Types</b>        | tissue homogenates or other biological fluids. |
| <b>Uniprot ID</b>          | <b>O60610</b>                                  |
| <b>Sensitivity</b>         | 0.051ng/mL                                     |
| <b>Entrez</b>              | <b>1729</b>                                    |
| <b>Expiration Date</b>     | Please enquire.                                |