

Product Datasheet

Human DEFa3 ELISA Kit (orb1199214)

Description	Human Defensin Alpha 3, Neutrophil Specific (DEFa3) ELISA Kit
Reactivity	Human
Range	46.8-3000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P59666
Sensitivity	16.47pg/mL
Entrez	1668
Expiration Date	Please enquire.