

Product Datasheet

Human DDAH2 ELISA Kit (orb1199219)

Description	Human Dimethylarginine Dimethylaminohydrolase 2 (DDAH2) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O95865
Sensitivity	0.119ng/mL
Entrez	23564
Expiration Date	Please enquire.