

Product Datasheet

Human DBI ELISA Kit (orb1199227)

Description	Human Diazepam Binding Inhibitor (DBI) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P07108
Sensitivity	0.054ng/mL
Entrez	1622
Expiration Date	Please enquire.