

Product Datasheet

Human CYP1A2 ELISA Kit (orb1199245)

Description	Human Cytochrome P450 1A2 (CYP1A2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P05177
Sensitivity	0.051ng/mL
Entrez	1544
Expiration Date	Please enquire.