

## Product Datasheet

### Human CYP11B1 ELISA Kit (orb1199247)

|                            |                                                                  |
|----------------------------|------------------------------------------------------------------|
| <b>Description</b>         | Human Cytochrome P450 11B1<br>(CYP11B1) ELISA Kit                |
| <b>Reactivity</b>          | Human                                                            |
| <b>Range</b>               | 0.156-10ng/mL                                                    |
| <b>Tested Applications</b> | ELISA                                                            |
| <b>Note</b>                | For research use only                                            |
| <b>Application notes</b>   | sandwich                                                         |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates or<br>other biological fluids. |
| <b>Uniprot ID</b>          | <b>P15538</b>                                                    |
| <b>Sensitivity</b>         | 0.074ng/mL                                                       |
| <b>Entrez</b>              | <b>1584</b>                                                      |
| <b>Expiration Date</b>     | Please enquire.                                                  |