

## Product Datasheet

### Human CYP11A1 ELISA Kit (orb1199248)

|                            |                                                   |
|----------------------------|---------------------------------------------------|
| <b>Description</b>         | Human Cytochrome P450 11A1<br>(CYP11A1) ELISA Kit |
| <b>Reactivity</b>          | Human                                             |
| <b>Range</b>               | 0.312-20ng/mL                                     |
| <b>Tested Applications</b> | ELISA                                             |
| <b>Note</b>                | For research use only                             |
| <b>Application notes</b>   | sandwich                                          |
| <b>Sample Types</b>        | tissue homogenates or other<br>biological fluids. |
| <b>Uniprot ID</b>          | <b>P05108</b>                                     |
| <b>Sensitivity</b>         | 0.133ng/mL                                        |
| <b>Entrez</b>              | <b>1583</b>                                       |
| <b>Expiration Date</b>     | Please enquire.                                   |