

Product Datasheet

Human CSN1 ELISA Kit (orb1199274)

Description	Human Casein Alpha (CSN1) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates, breast milk or other biological fluids.
Uniprot ID	P47710
Sensitivity	0.55ng/mL
Entrez	1446
Expiration Date	Please enquire.