

Product Datasheet

Human CRIM1 ELISA Kit (orb1199286)

Description	Human Cysteine Rich Transmembrane BMP Regulator 1 (CRIM1) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q9NZV1
Sensitivity	0.133ng/mL
Entrez	51232
Expiration Date	Please enquire.