

Product Datasheet

Human COPS2 ELISA Kit (orb1199318)

Description	Human COP9 Constitutive Photomorphogenic Homolog Subunit 2 (COPS2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P61201
Sensitivity	0.065ng/mL
Entrez	9318
Expiration Date	Please enquire.