

Product Datasheet

Human CNTNAP2 ELISA Kit (orb1199336)

Description	Human Contactin Associated Protein Like Protein 2 (CNTNAP2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9UHC6
Sensitivity	0.048ng/mL
Entrez	26047
Expiration Date	Please enquire.