

Product Datasheet

Human CDH17 ELISA Kit (orb1199386)

Description	Human Cadherin 17 (CDH17) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q12864
Sensitivity	0.054ng/mL
Entrez	1015
Expiration Date	Please enquire.