

Product Datasheet

Human CCbL1 ELISA Kit (orb1199395)

Description	Human Cysteine Conjugate Beta Lyase, Cytoplasmic (CCbL1) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q16773
Sensitivity	27pg/mL
Entrez	883
Expiration Date	Please enquire.