

Product Datasheet

Human BSDL ELISA Kit (orb1199432)

Description	Human Lipase, Bile Salt Dependent (BSDL) ELISA Kit
Reactivity	Human
Range	12.5-800ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	X6R868
Sensitivity	5.1ng/mL
Entrez	1056
Expiration Date	Please enquire.