

Product Datasheet

Human BECN1 ELISA Kit (orb1199446)

Description	Human Beclin 1 (BECN1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q14457
Sensitivity	0.059ng/mL
Entrez	8678
Expiration Date	Please enquire.