

Product Datasheet

Human BAMBI ELISA Kit (orb1199455)

Description	Human BMP And Activin Membrane Bound Inhibitor Homolog (BAMBI) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q13145
Sensitivity	0.064ng/mL
Entrez	25805
Expiration Date	Please enquire.