

Product Datasheet

Human BAG3 ELISA Kit (orb1199457)

Description	Human Bcl2 Associated Athanogene 3 (BAG3) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	O95817
Sensitivity	0.069ng/mL
Entrez	9531
Expiration Date	Please enquire.