

Product Datasheet

Human ATP4b ELISA Kit (orb1199464)

Description	Human ATPase, H ⁺ /K ⁺ Exchanging Beta Polypeptide (ATP4b) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P51164
Sensitivity	0.059ng/mL
Entrez	496
Expiration Date	Please enquire.