

Product Datasheet

Human ATP2C1 ELISA Kit (orb1199465)

Description	Human Calcium-Transporting ATPase Type 2C Member 1 (ATP2C1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P98194
Sensitivity	0.077ng/mL
Entrez	27032
Expiration Date	Please enquire.