

Product Datasheet

Human ATP1a1 ELISA Kit (orb1199470)

Description	Human ATPase, Na ⁺ /K ⁺ Transporting Alpha 1 Polypeptide (ATP1a1) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P05023
Sensitivity	26.6pg/mL
Entrez	476
Expiration Date	Please enquire.