

Product Datasheet

Human APOL1 ELISA Kit (orb1199504)

Description	Human Apolipoprotein L (APOL1) ELISA Kit
Reactivity	Human
Range	0.468-30ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	014791
Sensitivity	0.19ng/mL
Entrez	5842
Expiration Date	Please enquire.