

Product Datasheet

Human AOA^H ELISA Kit (orb1199514)

Description	Human Acyloxyacyl Hydrolase (AOAH) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P28039
Sensitivity	0.059ng/mL
Entrez	313
Expiration Date	Please enquire.